



Isocitrate Microplate Assay Kit

Catalog # AS0157

Detection and Quantification of Total Isocitrate Content in Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids Samples.

This instruction must be read in its entirety before using this product.

For research use only, Not for use in diagnostic procedures.

Contact information:

Tel:+1 (301) 446-2499 Fax:+1 (301) 446-2413

Email:techcn@signalwayantibody.com Web:www.sabbiotech.com

I. INTRODUCTION.....	2
II. KIT COMPONENTS.....	3
III. MATERIALS REQUIRED BUT NOT PROVIDED.....	3
VI. SAMPLE PREPARATION.....	4
V. ASSAY PROCEDURE.....	5
VI. CALCULATION.....	6
VII. TYPICAL DATA.....	7
VIII. TECHNICAL SUPPORT.....	7
IX. NOTES.....	7

I. INTRODUCTION

Isocitrate is a substrate in the citric acid (TCA) cycle. Isocitrate is formed by the isomerization of citrate catalyzed by the enzyme aconitase. Isocitrate is oxidized by isocitrate dehydrogenase producing α -ketoglutarate and generating NADPH. Isocitrate is commonly found in many fruits and vegetables and their processed products. Industrially, isocitrate is used as a marker to identify the quality and purity of fruit juices.

Isocitrate Microplate Assay Kit provides a simple and direct procedure for measuring isocitrate content in a variety of samples. The NADPH converts the dye to an intense violet color with an absorption maximum at 570 nm. The increase in absorbance at 570 nm is directly proportional to the isocitrate concentration.

II. KIT COMPONENTS

Component	Volume	Storage
96-Well Microplate	1 plate	
Assay Buffer	30ml x 4	4 °C
Enzyme Diluent	10 ml x 1	4 °C
Enzyme	Powder x 1	-20 °C
Dye Reagent	Powder x 1	4 °C
Dye Reagent Diluent	10 ml x 1	4 °C
Standard	Powder x 1	-20 °C
Technical Manual	1 Manual	

Note:

Enzyme: add 8 ml Enzyme Diluent to dissolve before use.

Standard: add 1 ml distilled water to dissolve before use; then add 0.25 ml into 0.75 ml distilled water, the concentration will be 500 $\mu\text{mol/L}$.

III. MATERIALS REQUIRED BUT NOT PROVIDED

1. Microplate reader to read absorbance at 570 nm
2. Distilled water
3. Pipettor
4. Pipette tips
5. Mortar
6. Centrifuge
7. Timer
8. Ice

IV. SAMPLE PREPARATION

1. For cell and bacteria samples

Collect cell or bacteria into centrifuge tube, discard the supernatant after centrifugation, add 1 ml Assay buffer for 5×10^6 cell or bacteria, sonicate (with power 20%, sonication 3s, interval 10s, repeat 30 times); centrifuged at 8000g 4°C for 10 minutes, take the supernatant into a new centrifuge tube and keep it on ice for detection.

2. For tissue samples

Weigh out 0.1 g tissue, homogenize with 1 ml Assay buffer on ice, centrifuged at 8000g 4°C for 10 minutes, take the supernatant into a new centrifuge tube and keep it on ice for detection.

3. For liquid samples

Detect directly.

V. ASSAY PROCEDURE

Warm all reagents to room temperature before use.

Add following reagents into the microplate:

Reagent	Sample	Standard	Blank
Sample	20 μ l	--	--
Enzyme	80 μ l	--	--
Incubate at room temperature for 10 minutes.			
Standard	--	100 μ l	--
Distilled water	--	--	100 μ l
Dye Reagent	100 μ l	100 μ l	100 μ l
Mix, wait for 2 minutes, measured at 570 nm and record the absorbance.			

VI. CALCULATION

1. According to the protein concentration of sample

$$\text{Isocitrate } (\mu\text{mol/mg}) = \frac{(C_{\text{Standard}} \times V_{\text{Standard}}) \times (OD_{\text{Sample}} - OD_{\text{Blank}})}{(OD_{\text{Standard}} - OD_{\text{Blank}}) \times (V_{\text{Sample}} \times C_{\text{Protein}})}$$

$$= 2.5 \times (OD_{\text{Sample}} - OD_{\text{Blank}}) / (OD_{\text{Standard}} - OD_{\text{Blank}}) / C_{\text{Protein}}$$

2. According to the weight of sample

$$\text{Isocitrate } (\mu\text{mol/g}) = \frac{(C_{\text{Standard}} \times V_{\text{Standard}}) \times (OD_{\text{Sample}} - OD_{\text{Blank}})}{(OD_{\text{Standard}} - OD_{\text{Blank}}) \times (V_{\text{Sample}} \times W / V_{\text{Assay}})}$$

$$= 2.5 \times (OD_{\text{Sample}} - OD_{\text{Blank}}) / (OD_{\text{Standard}} - OD_{\text{Blank}}) / W$$

3. According to the volume of serum or plasma

$$\text{Isocitrate } (\mu\text{mol/ml}) = \frac{(C_{\text{Standard}} \times V_{\text{Standard}}) \times (OD_{\text{Sample}} - OD_{\text{Blank}})}{(OD_{\text{Standard}} - OD_{\text{Blank}}) \times V_{\text{Sample}}}$$

$$= 2.5 \times (OD_{\text{Sample}(130\text{S})} - OD_{\text{Sample}(10\text{S})}) / (OD_{\text{Standard}} - OD_{\text{Blank}})$$

C_{Standard} : the standard concentration, 500 $\mu\text{mol/L}$ = 0.5 $\mu\text{mol/ml}$;

V_{Standard} : the volume of standard, 100 μl = 0.1 ml;

C_{Protein} : the protein concentration, mg/ml;


W : the weight of sample, g;

V_{Sample} : the volume of sample, 20 μl = 0.02 ml;

V_{Assay} : the volume of Assay buffer, 1 ml.

VII. TYPICAL DATA

The standard curve is for demonstration only. A standard curve must be run with each assay.



Detection Range: 5µmol/L - 500 µmol/L

VIII. TECHNICAL SUPPORT

For troubleshooting, information or assistance, please go online to www.sabbiotech.cn or contact us at techcn@signalwayantibody.com

IX. NOTES