

PASK antibody

Catalog No: #22236

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	PASK antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 930 and 1311 of Human PASK
Target Name	PASK
Accession No.	NCBI Gene ID: 23178NCBI mRNA#: BC063585NCBI Protein#: AAH63585
Concentration	0.9mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 143kd

Western blotting: 1:500-1:3000

Images



Sample (30ug whole cell lysate)
A: Raji
5% SDS PAGE
Primary antibody diluted at 1: 1000

Background

PAS domains regulate the function of many intracellular signaling pathways in response to both extrinsic and intrinsic stimuli. PASK is an evolutionarily conserved protein present in yeast, flies, and mammals.[supplied by OMIM]

Note: This product is for in vitro research use only and is not intended for use in humans or animals.