

RFP Mouse Polyclonal Antibody

Catalog No: #T511

Package Size: #T511-1 50ul #T511-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	RFP Mouse Polyclonal Antibody
Host Species	Mouse
Clonality	Polyclonal
Purification	Affinity purification using immunogen.
Applications	WB IHC IP
Specificity	The RFP tag antibody can recognize RFP tag fusion proteins.
Target Name	RFP
Concentration	1.0mg/ml
Formulation	Mouse IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:2000~1:5000

Images



Western blot analysis of RFP Recombinant protein, using #T511 diluted at 1) 1:5,000, 2) 1:10,000.

Background

The red fluorescent protein cloned from *Discosoma coral* (DsRed or drFP583) holds great promise for biotechnology and cell biology as a spectrally distinct companion or substitute for the green fluorescent protein (GFP) from the *Aequorea jellyfish*.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.