

## IL13RA2 antibody

Catalog No: #38341

Package Size: #38341-1 50ul #38341-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	IL13RA2 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were purified by affinity purification using immunogen.
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total IL13RA2 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant protein of human IL13RA2.
Target Name	IL13RA2
Other Names	CT19; IL-13R; IL13BP; CD213A2;
Accession No.	Swiss-Prot#: Q14627NCBI Gene ID: 3598
SDS-PAGE MW	44kd
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

WB 1:500 - 1:2000 IHC 1:50 - 1:200

## Images



Western blot analysis of extracts of various cell lines, using IL13RA2 antibody at 1:1000 dilution.

## Background

The protein encoded by this gene is closely related to IL13RA1, a subunit of the interleukin 13 receptor complex. This protein binds IL13 with high affinity, but lacks cytoplasmic domain, and does not appear to function as a signal mediator. It is reported to play a role in the internalization of IL13.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.