ZBTB38 antibody

Catalog No: #38964

Package Size: #38964-1 50ul #38964-2 100ul



Orders: order@signalwayantibody.com Support: tech@signalwayantibody.com


D	cription
1 100	crintion
	UIDUUUI

Becomption	
Product Name	ZBTB38 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were purified by affinity purification using immunogen.
Applications	WB
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total ZBTB38 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant protein of human ZBTB38.
Target Name	ZBTB38
Other Names	CIBZ; ZNF921; PPP1R171;
Accession No.	Swiss-Prot#: Q8NAP3NCBI Gene ID: 253461
SDS-PAGE MW	134kd
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02%
	sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using ZBTB38 antibody at 1:1000 dilution.

Background

The protein encoded by this gene is a zinc finger transcriptional activator that binds methylated DNA. The encoded protein can form homodimers or heterodimers through the zinc finger domains. In mouse, inhibition of this protein has been associated with apoptosis in some cell types.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.