

Recombinant Human Thymidine Kinase,
Cytosolic,TK1 (C-6His)

Catalog No: #AP76035



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Thymidine Kinase, Cytosolic,TK1 (C-6His)
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P04183
Calculated MW	26.5kDa
Target Sequence	MSCINLPTVLPGPSKTRGQIQVILGPMFSGKSTELMRRVRRFQIAQYKCLVIKYAKDTRYSSSFCTHDRNTME ALPACLLRDVAQEALGVAVIGIDEGQFFPDIVEFCEAMANAGKTVIVAALDGTFRKPFGAILNLVPLAESVVKL TAVCMECFREAAAYTKRLGTEKEVEVIGGADKYHVCRLCYFKKASGQPAGPDNKENCVPVPGKPGEAVAARKL FAPQQILQCSPANVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of PBS,pH7.4.
Storage	Store at < -70°C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.