

# EBFP Monoclonal Antibody

Catalog No: #T615



Package Size: #T615-1 50ul #T615-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	EBFP Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	The antibody was affinity-purified from antiserum by affinity-chromatography
Applications	WB
Specificity	The antibody detects EBFP & EBFP Tag protein.
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Western blotting: 1:5000

## Images



Western blot analysis of recombinant EBFP protein, diluted at 1:5000.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.