

UQCC2 Antibody

Catalog No: #47272

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	UQCC2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu, Ms, Rt
Specificity	The antibody detects endogenous levels of total UQCC2 protein.
Immunogen Type	Peptide
Immunogen Description	Fusion protein of human UQCC2
Target Name	UQCC2
Other Names	M19; Cbp6; MNF1; C6orf125; C6orf126; bA6B20.2
Accession No.	Swiss-Prot#:Q9BRT2NCBI Gene ID:84300Gene Accssion:BC006007
Concentration	1.1
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

IHC dilution:1: 40-200

Images



The image is immunohistochemistry of paraffin-embedded Human breast cancer tissue using 47272(UQCC2 Antibody) at dilution 1/60.(Original magnification: 200)

Background

This gene encodes a nucleoid protein localized to the mitochondria inner membrane. The encoded protein affects regulation of insulin secretion, mitochondrial ATP production, and myogenesis through modulation of mitochondrial respiratory chain activity.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.