

TGFβ3 Polyclonal Antibody

Catalog No: #41496

Package Size: #41496-1 50ul #41496-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TGFβ3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA IF
Species Reactivity	Hu Ms Rt Mk
Specificity	TGFβ3 Polyclonal Antibody detects endogenous levels of TGFβ3 protein.
Immunogen Description	Synthesized peptide derived from the Internal region of human TGFβ3.
Target Name	TGFβ3
Other Names	TGFB3; Transforming growth factor beta-3; TGF-beta-3
Accession No.	Swiss-Prot: P10600NCBI Gene ID: 7043
SDS-PAGE MW	13kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/10000.

Not yet tested in other applications.

IF 1:50-200

Images



Western Blot analysis of COS7 cells using TGFβ3 Polyclonal Antibody

Published Papers

et al., A Dermal Equivalent Engineered with TGFβ1 and β2 Expressing Bone Marrow Stromal Cells and Amniotic Membrane: Cosmetic Healing of Full-thickness Skin Wounds in Rats. In Oncotarget on 2016 Feb 2 by Pei Chen, Na Yu, et al..PMID:26756218, (2016)
[PMID:26756218](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.