

Rabbit IgG

Catalog No: #IG401



Package Size: #IG401-1 100mg #IG401-2 200mg #IG401-3 1g

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|---------------|---|
| Product Name | Rabbit IgG |
| Purification | Affinity chromatography or caprylic acid-ammonium sulfate,β 95% |
| Conjugates | None |
| Concentration | mg/ml |
| Storage | -20°C |

Product Description

Source: Normal, healthy, unimmunized rabbit serum

Appearance: Liquid

Filtration: 0.2um

Buffer :10mM Phosphate buffered salinePBSpH7.4

Preservative: None

Stability: Stable for 24 months

Precautions: This product can be kept at 2-8/ C for maximum 30 days. Avoid repeated freezing and thawing. If multiple usages are needed, it is recommended to be stored separately in smaller volumes.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.