

Chicken IgY

Catalog No: #IG601

Package Size: #IG601-1 100mg #IG601-2 200mg #IG601-3 1g

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Chicken IgY
Purification	Water dilution- sodium sulfate ,β 85%
Conjugates	None
Concentration	β 20 mg/ml
Storage	-20°C

Product Description

Source: Normal chicken egg yolk

Appearance: Liquid

Filtration: 0.2um

Buffer :10mM Phosphate buffered salinePBSpH7.4

Preservative: None

Stability: Stable for 24 months

Precautions: This product can be kept at 2-8/ C for maximum 30 days. Avoid repeated freezing and thawing. If multiple usages are needed, it is recommended to be stored separately in smaller volumes.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.