

Recombinant Lassa virus Nucleoprotein(N)

Catalog No: #AP77647

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Lassa virus Nucleoprotein(N)
Brief Description	Recombinant Protein
Host Species	Yeast
Purification	Greater than 85% as determined by SDS-PAGE.
Immunogen Description	Expression Region:1-570aaSequence Info:Full Length
Target Name	N
Other Names	Nucleocapsid protein Protein N
Accession No.	P04935
Calculated MW	62.9 kDa
Tag Info	NO-tagged
Target Sequence	MSASKEVRSLWTQSLRRELSGYCSNIKLQVVKDAQALLHGLDFSEVSNVQRLMRKQKRDDGDLKRLRDLN QAVNNLVELKSTQQKSVLRVGTLSSDLLILAADLEKLKSKVTRTERPLSSGVYMGNLSSQQLDQRRALLNMI GMTGVSGGGKGASNGIVRVWDVKNAELLNNQFGTMPSLTLACLTQGGQVDLNDVQALTDLGLIYTAKYPNS SDLDRLSQSHPIILNMIDTKKSSLNISGYNFSLGAAVKAGACMLDGGNMLETIKVSPQTMKGILKSILKVKKSLGM FVSDTPGERNPYENILYKICLSGDGWSYIASRTSIVGRAWENTVVDLEQDNKPQKIGNGGSNKSLSAGFAAG LIYSQLMTLKDFKCFNLIPNAKTWMDIEGRPEDPVEIALYQPSSGCYVHFFREPTDLKQFKQDAKYSHGRDVT DLFAAQPGLTSAVIEALPRNMVITCQGSSEDIRKLLESQGRDIDITLSKADSRKFENAVVDQFKDLCHMHT GVVVEKKKRGKKEITPHCALMDCIMFDAAVSGGLDAKVLRVVLPMDVFRTPKVVVL
Formulation	Tris-based buffer50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Application Details

tissue_specificity:Homo sapiens (Human)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.