

# Guanine nucleotide-binding protein G(s) subunit alpha isoforms short

Catalog No: #AP79721

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Guanine nucleotide-binding protein G(s) subunit alpha isoforms short
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Alternative gene product encoded by XL-exon,GNAS1
Accession No.	Uniprot ID:P84996
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.