

Bifunctional heparan sulfate
N-deacetylase/N-sulfotransferase 2

Catalog No: #AP79729



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Bifunctional heparan sulfate N-deacetylase/N-sulfotransferase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Glucosaminyl N-deacetylase/N-sulfotransferase 2,NDST-2,N-heparan sulfate sulfotransferase 2,N-HSST 2,HSST2
Accession No.	Uniprot ID:P52849
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.