

Proteasome subunit beta type-9

Catalog No: #AP79731



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Proteasome subunit beta type-9
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Low molecular mass protein 2,Macropain chain 7,Multicatalytic endopeptidase complex chain 7,Proteasome chain 7,Proteasome subunit beta-1i,Really interesting new gene 12 protein,Lmp2,Ring12
Accession No.	Uniprot ID:P28077
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.