

Paired immunoglobulin-like type 2 receptor beta

Catalog No: #AP79733



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Description

Product Name	Paired immunoglobulin-like type 2 receptor beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Activating receptor PILR-beta,Cell surface receptor FDFACT,FDFACT,PP1551
Accession No.	Uniprot ID:Q9UKJ0
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.