

Heat shock protein beta-1

Catalog No: #AP79738



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Heat shock protein beta-1 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Hu |
| Other Names | HspB1,28 kDa heat shock protein,Estrogen-regulated 24 kDa protein,Heat shock 27 kDa protein,HSP 27,Stress-responsive protein 27,SRP27,HSP27,HSP28 |
| Accession No. | Uniprot ID:P04792 |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Note: This product is for in vitro research use only and is not intended for use in humans or animals.