

Jouberin

Catalog No: #AP80436

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Jouberin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Abelson helper integration site 1 protein,AHI-1
Accession No.	Uniprot ID:Q8K3E5
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.