

# Acetyl-Histone H3 (K9) Monoclonal Antibody(2E7)

Catalog No: #HW205

Package Size: #HW205-1 50ul #HW205-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	Acetyl-Histone H3 (K9) Monoclonal Antibody(2E7)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	2.00E+07
Isotype	IgG1
Applications	IHC
Species Reactivity	Hu Rt Ms
Specificity	Acetyl-Histone H3 (K9) Monoclonal Antibody(2E7) detects endogenous levels of HIST1H3A
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide
Target Name	Histone H3
Other Names	HIST1H3A
Calculated MW	15kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


## Application Details

Immunohistochemistry: 1:100-1:200

## Images



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma Tissue using Acetyl Histone H3 K9 Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Human Stomach Tissue using Acetyl Histone H3 K9 Mouse mAb diluted at 1:200.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.