

DYRK2/4 (Phospho-Tyr382/264) Antibody

Catalog No: #12498

Package Size: #12498-1 50ul #12498-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|---|
| Product Name | DYRK2/4 (Phospho-Tyr382/264) Antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide. |
| Applications | WB |
| Species Reactivity | Hu Ms Rt |
| Specificity | DYRK2/4 (Phospho-Tyr382/264) Antibody detects endogenous levels of DYRK2/4 only when phosphorylated at Tyr382/264 |
| Immunogen Type | Peptide |
| Immunogen Description | A synthesized peptide derived from human DYRK2/4 (Phospho-Tyr382/264) |
| Target Name | DYRK2/4 |
| Modification | Phospho |
| Other Names | DYRK2 |
| Accession No. | Swiss-Prot#: Q92630/Q9NR20NCBI Gene ID: 8445/8798 |
| Target Species | human |
| Calculated MW | 60kd |
| Concentration | 1.0mg/ml |
| Formulation | Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Storage | Store at -20°C |

Application Details

Western blotting: 1:1000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.