

WAP four-disulfide core domain protein 2

Catalog No: #AP89200



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	WAP four-disulfide core domain protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Epididymal secretory protein E4
Accession No.	Uniprot ID:Q8MI69
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.