

MCP-1 Antibody

Catalog No: #45071



Package Size: #45071-1 50ul #45071-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MCP-1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	MCP-1 Antibody detects endogenous levels of total MCP-1
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human MCP-1
Target Name	MCP-1
Other Names	CCL2, C-C motif chemokine 2, Chemokine (C-C motif) ligand 2, GDCF-2, HC11, HSMCR30, MCP1, SCYA2, Small-inducible cytokine A2, Monocyte secretory protein JE, MCAF, MCP-1, Monocyte chemotactic protein 1, SMC-CF
Accession No.	Swiss-Prot#: P13500NCBI Gene ID:6347
Calculated MW	14 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Published Papers

el at., Screening of anti ζ Ahronic nonbacterial prostatitis activity of different extractions of the aerial part of Glycyrrhiza uralensis, and network pharmacology research. In Biomed Rep on 2021 Dec by Haifan Liu, Jie Cui,et al..PMID:34667596, , (2021)

[PMID:34667596](https://pubmed.ncbi.nlm.nih.gov/34667596/)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.