

OR6N1 Antibody

Catalog No: #46112



Package Size: #46112-1 50ul #46112-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	OR6N1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Applications	WB
Species Reactivity	Hu
Specificity	OR6N1 Antibody detects endogenous levels of total OR6N1
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human OR6N1
Target Name	OR6N1
Other Names	OR6N1, Olfactory receptor OR1-22, Olfactory receptor 6N1, OR1-22
Accession No.	Swiss-Prot#: Q8NGY5NCBI Gene ID:128372
Calculated MW	35 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis OR6N1 using 293 whole cell lysates

Note: This product is for in vitro research use only and is not intended for use in humans or animals.