

# HTR1F Polyclonal Antibody

Catalog No: #28691



Package Size: #28691-1 50ul #28691-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	HTR1F Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Hu,Ms,Rt
Immunogen Description	Recombinant fusion protein of human HTR1F (NP_000857.1).
Other Names	HTR1F; 5-HT-1F; 5-HT1F; 5HT6; HTR1EL; MR77; 5-hydroxytryptamine receptor 1F
Accession No.	Swiss-Prot#:P30939NCBI Gene ID:3355
Calculated MW	42kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

## Application Details

WB 1:500 - 1:2000

## Images



Western blot analysis of extracts of various cell lines, using HTR1F antibody.

## Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.