

HAO2 Polyclonal Antibody

Catalog No: #28846

Package Size: #28846-1 50ul #28846-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	HAO2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Ms
Immunogen Description	Recombinant fusion protein of human HAO2 (NP_057611.1).
Other Names	HAO2; GIG16; HAOX2; hydroxyacid oxidase 2
Accession No.	Swiss-Prot#:Q9NYQ3NCBI Gene ID:51179
Calculated MW	38kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of Mouse kidney, using HAO2 antibody.

Background

This gene is one of three related genes that have 2-hydroxyacid oxidase activity. The encoded protein localizes to the peroxisome has the highest activity toward the substrate 2-hydroxypalmitate. Alternative splicing results in multiple transcript variants.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.