

OTUD7A Polyclonal Antibody

Catalog No: #29345



Package Size: #29345-1 50ul #29345-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	OTUD7A Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Hu
Immunogen Description	A synthetic peptide of mouse OTUD7A (NP_570950.1).
Other Names	OTUD7A; C15orf16; C16ORF15; CEZANNE2; OTUD7; OTU deubiquitinase 7A
Accession No.	Swiss-Prot#:_NCBI Gene ID:170711
Calculated MW	105kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:200 - 1:2000

Images



Western blot analysis of extracts of U-251MG cells, using OTUD7A antibody.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.