

HNF4G Polyclonal Antibody

Catalog No: #29544

Package Size: #29544-1 50ul #29544-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	HNF4G Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Hu,Ms,Rt
Immunogen Description	Recombinant fusion protein of human HNF4G (NP_004124.4).
Other Names	HNF4G; NR2A2; NR2A3; hepatocyte nuclear factor 4-gamma
Accession No.	Swiss-Prot#:Q14541NCBI Gene ID:3174
Calculated MW	50kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB 1:500 - 1:2000IF 1:50 - 1:100

Images



Western blot analysis of extracts of various cell lines, using HNF4G antibody.



Immunofluorescence analysis of U-2 OS cells using HNF4G Polyclonal Antibody.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.