

KIAA1429 Polyclonal Antibody

Catalog No: #29774



Package Size: #29774-1 50ul #29774-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	KIAA1429 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human KIAA1429 (NP_056311.2).
Other Names	KIAA1429; MSTP054; fSAP121
Accession No.	Swiss-Prot#:Q69YN4NCBI Gene ID:25962
Calculated MW	208kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB 1:500 - 1:2000IHC 1:50 - 1:100IF 1:50 - 1:100


Images




Western blot analysis of extracts of various cell lines, using KIAA1429 antibody at 1:1000 dilution.




Immunohistochemistry of paraffin-embedded rat brain using KIAA1429 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded human thyroid cancer using KIAA1429 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded human placenta using KIAA1429 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded mouse brain using KIAA1429 antibody at dilution of 1:100 (40x lens).

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.