

# OGFOD1 Polyclonal Antibody

Catalog No: #29873



Package Size: #29873-1 50ul #29873-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

|                       |                                                                   |
|-----------------------|-------------------------------------------------------------------|
| Product Name          | OGFOD1 Polyclonal Antibody                                        |
| Host Species          | Rabbit                                                            |
| Clonality             | Polyclonal                                                        |
| Isotype               | IgG                                                               |
| Purification          | Affinity purification                                             |
| Applications          | WB                                                                |
| Species Reactivity    | Hu                                                                |
| Immunogen Description | Recombinant fusion protein of human OGFOD1 (NP_060703.3).         |
| Other Names           | OGFOD1; TPA1; prolyl 3-hydroxylase OGFOD1                         |
| Accession No.         | Swiss-Prot#:Q8N543NCBI Gene ID:55239                              |
| Calculated MW         | 63kDa                                                             |
| Formulation           | Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4. |
| Storage               | Store at -20°C                                                    |

## Application Details

WB 1:500 - 1:2000

## Images



Western blot analysis of extracts of HeLa cells, using OGFOD1 antibody.

## Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.