

DHX57 Polyclonal Antibody

Catalog No: #29895



Package Size: #29895-1 50ul #29895-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	DHX57 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Hu,Rt
Immunogen Description	Recombinant fusion protein of human DHX57 (NP_945314.1).
Other Names	DHX57; DDX57; DExH-box helicase 57
Accession No.	Swiss-Prot#:Q6P158NCBI Gene ID:90957
Calculated MW	150-155kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using DHX57 antibody.

Background

Published Papers

el at., High Levels of DEAH-Box Helicases Relate to Poor Prognosis and Reduction of DHX9 Improves Radiosensitivity of Hepatocellular Carcinoma. In Front Oncol on 2022 Jun 22 by Xi Chen, Letao Lin, et al..PMID: 35814441, , (2022)
[PMID:35814441](https://pubmed.ncbi.nlm.nih.gov/35814441/)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.