

MRGPRF(Grp168)Polyclonal Antibody

Catalog No: #30338

Package Size: #30338-1 50ul #30338-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MRGPRF(Grp168)Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Hu
Immunogen Description	A synthetic peptide of human MRGPRFα(1-168)Grp168α(1-168).
Accession No.	Swiss-Prot#:Q96AM1NCBI Gene ID:116535
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.