

# OSBPL3 Polyclonal Antibody

Catalog No: #30591



Package Size: #30591-1 50ul #30591-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	OSBPL3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Hu,Ms,Rt
Immunogen Description	Recombinant fusion protein of human OSBPL3 (NP_056365.1).
Other Names	OSBPL3; ORP-3; ORP3; OSBP3; oxysterol binding protein like 3
Accession No.	Swiss-Prot#:Q9H4L5NCBI Gene ID:26031
Calculated MW	101kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

## Application Details

WB 1:200 - 1:2000IHC 1:50 - 1:200

## Images



Western blot analysis of extracts of various cell lines, using OSBPL3 antibody.

## Background

This gene encodes a member of the oxysterol-binding protein (OSBP) family, a group of intracellular lipid receptors. Most members contain an N-terminal pleckstrin homology domain and a highly conserved C-terminal OSBP-like sterol-binding domain. The encoded protein is involved in the regulation of cell adhesion and organization of the actin cytoskeleton. Alternative splicing results in multiple transcript variants.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.