

QRSL1 Polyclonal Antibody

Catalog No: #31573



Package Size: #31573-1 50ul #31573-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	QRSL1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Hu,Ms
Immunogen Description	Recombinant fusion protein of human QRSL1 (NP_060762.3).
Other Names	QRSL1; GatA; glutaminyl-tRNA synthase (glutamine-hydrolyzing)-like 1
Accession No.	Swiss-Prot#:Q9H0R6NCBI Gene ID:55278
Calculated MW	57kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using QRSL1 antibody.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.