

OSM Polyclonal Antibody

Catalog No: #46796



Package Size: #46796-1 50ul #46796-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	OSM Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous OSM
Immunogen Type	protein
Immunogen Description	Synthetic peptide from human protein at AA range: 50-100
Target Name	OSM
Other Names	Oncostatin-M (OSM)
Accession No.	Swiss-Prot:P13725NCBI Gene ID:5008
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C


Application Details

Immunohistochemistry: IHC 1:50-300


Images




Immunohistochemical analysis of paraffin-embedded human-brain, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-brain, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-colon, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-colon, antibody was diluted at 1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.