

OTUD6B Conjugated Antibody

Catalog No: #C28608

Package Size: #C28608-AF350 100ul #C28608-AF405 100ul #C28608-AF488 100ul #C28608-AF555 100ul #C28608-AF594 100ul #C28608-AF647 100ul #C28608-AF680 100ul #C28608-AF750 100ul #C28608-Biotin 100ul #C28608-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	OTUD6B Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB, IF
Species Reactivity	Hu,Ms,Rt
Immunogen Description	Recombinant fusion protein of human OTUD6B (NP_057107.3).
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	OTUD6B; CGI-77; DUBA-5; DUBA5; IDDFSDA; OTU domain containing 6B
Accession No.	Swiss-Prot#:Q8N6M0NCBI Gene ID:51633
Calculated MW	38kDa
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Background

This gene encodes a member of the ovarian tumor domain (OTU)-containing subfamily of deubiquitinating enzymes. Deubiquitinating enzymes are primarily involved in removing ubiquitin from proteins targeted for degradation. This protein may function as a negative regulator of the cell cycle in B cells.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.