

OGDHL Conjugated Antibody

Catalog No: #C28987

Package Size: #C28987-AF350 100ul #C28987-AF405 100ul #C28987-AF488 100ul #C28987-AF555 100ul #C28987-AF594 100ul
#C28987-AF647 100ul #C28987-AF680 100ul #C28987-AF750 100ul #C28987-Biotin 100ul #C28987-Conjugated 50ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	OGDHL Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB, IF
Species Reactivity	Hu,Ms,Rt
Immunogen Description	Recombinant fusion protein of human OGDHL (NP_060715.2).
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	OGDHL; oxoglutarate dehydrogenase-like
Accession No.	Swiss-Prot#:Q9ULD0NCBI Gene ID:55753
Calculated MW	120kDa
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Background

The protein encoded by this gene is similar to oxoglutarate dehydrogenase (OGDH) of the OGDH complex, which degrades glucose and glutamate. This gene encodes several isoforms, including some that appear to localize to mitochondria. The encoded protein down-regulates the AKT signaling cascade and can suppress the growth of cervical cancer cells.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.