

# OSTC Conjugated Antibody

Catalog No: #C29993

Package Size: #C29993-AF350 100ul #C29993-AF405 100ul #C29993-AF488 100ul #C29993-AF555 100ul #C29993-AF594 100ul #C29993-AF647 100ul #C29993-AF680 100ul #C29993-AF750 100ul #C29993-Biotin 100ul #C29993-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |                                                                                             |
|-----------------------|---------------------------------------------------------------------------------------------|
| Product Name          | OSTC Conjugated Antibody                                                                    |
| Host Species          | Rabbit                                                                                      |
| Clonality             | Polyclonal                                                                                  |
| Isotype               | IgG                                                                                         |
| Purification          | Affinity purification                                                                       |
| Applications          | WB, IF                                                                                      |
| Species Reactivity    | Hu,Ms,Rt                                                                                    |
| Immunogen Description | A synthetic synthetic peptide of human OSTC (NP_067050.1).                                  |
| Conjugates            | Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750                                      |
| Other Names           | OSTC; DC2; oligosaccharyltransferase complex subunit OSTC                                   |
| Accession No.         | Swiss-Prot#:Q9NRP0NCBI Gene ID:58505                                                        |
| Calculated MW         | Refer to figures                                                                            |
| Formulation           | 0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide |
| Storage               | Store at 4°C in dark for 6 months                                                           |

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.