

KIAA0513 antibody

Catalog No: #31860



Package Size: #31860-1 50ul #31860-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	KIAA0513 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu, Ms
Immunogen Description	Fusion protein of human KIAA0513
Target Name	KIAA0513
Accession No.	Swiss-Prot#: O60879NCBI Protein#: BC030280
Concentration	1.9mg/ml
Formulation	pH7.4 PBS, 0.05% NaN3, 40% Glycerol
Storage	Store at -20°C/1 year


Application Details

IHC dilution:1:50-1:300

Images



The image is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using (KIAA0513 Antibody) at dilution 1/50.



The image is immunohistochemistry of paraffin-embedded Human tonsil tissue using (KIAA0513 Antibody) at dilution 1/50.

Background

The function of this protein remains unknown.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.