

GRB14 Conjugated Antibody

Catalog No: #C31807



Package Size: #C31807-AF350 100ul #C31807-AF405 100ul #C31807-AF488 100ul #C31807-AF555 100ul #C31807-AF594 100ul #C31807-AF647 100ul #C31807-AF680 100ul #C31807-AF750 100ul #C31807-Biotin 100ul #C31807-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	GRB14 Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	WB, IF
Species Reactivity	Hu, Ms, Rt
Immunogen Description	Fusion protein of human GRB14
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Target Name	GRB14
Accession No.	Swiss-Prot#: Q9BTE6NCBI Protein#: BC053559
Calculated MW	61 kDa
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at -20°C/1 year

Application Details

WB: 1:50-1:200

IF: 1:50-1:200

Background

The product of this gene belongs to a small family of adapter proteins that are known to interact with a number of receptor tyrosine kinases and signaling molecules. This gene encodes a growth factor receptor-binding protein that interacts with insulin receptors and insulin-like growth-factor receptors. This protein likely has an inhibitory effect on receptor tyrosine kinase signaling and, in particular, on insulin receptor signaling.?

Note: This product is for in vitro research use only and is not intended for use in humans or animals.