

SURF6 Conjugated Antibody

Catalog No: #C31816



Package Size: #C31816-AF350 100ul #C31816-AF405 100ul #C31816-AF488 100ul #C31816-AF555 100ul #C31816-AF594 100ul #C31816-AF647 100ul #C31816-AF680 100ul #C31816-AF750 100ul #C31816-Biotin 100ul #C31816-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	SURF6 Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IF
Species Reactivity	Hu
Immunogen Description	Fusion protein of human SURF6
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Target Name	SURF6
Other Names	RRP14
Accession No.	Swiss-Prot#: Q68DC2NCBI Protein#: BC014878
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at -20°C/1 year

Application Details

IF:1:50-1:200

Background

This gene encodes a conserved protein that is localized to the nucleolus. The encoded protein may function as a nucleolar-matrix protein with nucleic acid-binding properties. There is a pseudogene for this gene on chromosome Y. Alternative splicing results in multiple transcript variants.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.