

ZNF266 Conjugated Antibody

Catalog No: #C31822

Package Size: #C31822-AF350 100ul #C31822-AF405 100ul #C31822-AF488 100ul #C31822-AF555 100ul #C31822-AF594 100ul #C31822-AF647 100ul #C31822-AF680 100ul #C31822-AF750 100ul #C31822-Biotin 100ul #C31822-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	ZNF266 Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IF
Species Reactivity	Hu
Immunogen Description	Fusion protein of human ZNF266
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Target Name	ZNF266
Other Names	HZF1
Accession No.	Swiss-Prot#: Q8IUR7NCBI Protein#: BC017407
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at -20°C/1 year

Application Details

IF:1:50-1:200

Background

This gene encodes a protein containing many tandem zinc-finger motifs. Zinc fingers are protein or nucleic acid-binding domains, and may be involved in a variety of functions, including regulation of transcription. This gene is located in a cluster of similar genes encoding zinc finger proteins on chromosome 19. Alternative splicing results in multiple transcript variants for this gene.?

Note: This product is for in vitro research use only and is not intended for use in humans or animals.