

ZNF300 Conjugated Antibody

Catalog No: #C31824

Package Size: #C31824-AF350 100ul #C31824-AF405 100ul #C31824-AF488 100ul #C31824-AF555 100ul #C31824-AF594 100ul
#C31824-AF647 100ul #C31824-AF680 100ul #C31824-AF750 100ul #C31824-Biotin 100ul #C31824-Conjugated 50ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ZNF300 Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IF
Species Reactivity	Hu
Immunogen Description	Fusion protein of human ZNF300
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Target Name	ZNF300
Accession No.	Swiss-Prot#: O15143NCBI Protein#: BC117248
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at -20°C/1 year

Application Details

IF: 1:50-1:200

Background

The protein encoded by this gene is a C2H2-type zinc finger DNA binding protein and likely transcriptional regulator. The function of this protein is not yet known. Three transcript variants encoding different isoforms have been found for this gene.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.