

ZNF410 Conjugated Antibody

Catalog No: #C31826

Package Size: #C31826-AF350 100ul #C31826-AF405 100ul #C31826-AF488 100ul #C31826-AF555 100ul #C31826-AF594 100ul #C31826-AF647 100ul #C31826-AF680 100ul #C31826-AF750 100ul #C31826-Biotin 100ul #C31826-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	ZNF410 Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IF
Species Reactivity	Hu
Immunogen Description	Fusion protein of human ZNF410
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Target Name	ZNF410
Other Names	APA1; APA-1
Accession No.	Swiss-Prot#: O75129NCBI Protein#: BC050683
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at -20°C/1 year

Application Details

IF:1:50-1:200

Background

ZNF410 is the transcription factor that activates transcription of matrix-remodeling genes such as MMP1 during fibroblast senescence.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.