

VWA2 Conjugated Antibody

Catalog No: #C31831

Package Size: #C31831-AF350 100ul #C31831-AF405 100ul #C31831-AF488 100ul #C31831-AF555 100ul #C31831-AF594 100ul #C31831-AF647 100ul #C31831-AF680 100ul #C31831-AF750 100ul #C31831-Biotin 100ul #C31831-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	VWA2 Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IF
Species Reactivity	Hu
Immunogen Description	Fusion protein of human VWA2
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Target Name	VWA2
Other Names	AMACO; NET42; CCSP-2
Accession No.	Swiss-Prot#: Q9UIR0NCBI Protein#: BC128588
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at -20°C/1 year

Application Details

IF:1:50-1:200

Background

This gene encodes a member of the von Willebrand factor A-like domain protein superfamily. The encoded protein is localized to the extracellular matrix and may serve as a structural component in basement membranes or in anchoring structures on scaffolds of collagen VII or fibrillin. This gene has been linked to type 1A diabetes and is a candidate serological marker for colon cancer.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.