

USP35 Conjugated Antibody

Catalog No: #C31832

Package Size: #C31832-AF350 100ul #C31832-AF405 100ul #C31832-AF488 100ul #C31832-AF555 100ul #C31832-AF594 100ul #C31832-AF647 100ul #C31832-AF680 100ul #C31832-AF750 100ul #C31832-Biotin 100ul #C31832-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	USP35 Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IF
Species Reactivity	Hu
Immunogen Description	Fusion protein of human USP35
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Target Name	USP35
Accession No.	Swiss-Prot#: Q9H0W9NCBI Protein#: BC131489
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at -20°C/1 year

Application Details

IF: 1:50-1:200

Background

USP35 (ubiquitin specific peptidase 35) is a 1,017 amino acid protein that belongs to the peptidase C19 family. Expressed in pancreas, testis and skeletal muscle, USP35 functions to catalyze the conversion of a ubiquitin C-terminal thioester and water to free ubiquitin and a thiol, a reaction that may play a role in signaling events throughout the cell.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.