

GATD3B Conjugated Antibody

Catalog No: #C31836

Package Size: #C31836-AF350 100ul #C31836-AF405 100ul #C31836-AF488 100ul #C31836-AF555 100ul #C31836-AF594 100ul #C31836-AF647 100ul #C31836-AF680 100ul #C31836-AF750 100ul #C31836-Biotin 100ul #C31836-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	GATD3B Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	WB, IF
Species Reactivity	Hu
Immunogen Description	Fusion protein of human GATD3B
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Target Name	GATD3B
Other Names	KNP-I
Accession No.	Swiss-Prot#: Q8N163 NCBI Protein#: BC003587
Calculated MW	28 kDa
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at -20°C/1 year

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Background

May play a role in promoting mitochondrial enlargement in cone photoreceptor cells.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.