

C3orf38 Conjugated Antibody

Catalog No: #C31837



Package Size: #C31837-AF350 100ul #C31837-AF405 100ul #C31837-AF488 100ul #C31837-AF555 100ul #C31837-AF594 100ul #C31837-AF647 100ul #C31837-AF680 100ul #C31837-AF750 100ul #C31837-Biotin 100ul #C31837-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	C3orf38 Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IF
Species Reactivity	Hu, Ms, Rt
Immunogen Description	Fusion protein of human C3orf38
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Target Name	C3orf38
Accession No.	Swiss-Prot#: Q9UPN4 NCBI Protein#: BC024188
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at -20°C/1 year

Application Details

IF: 1:50-1:200

Background

May be involved in apoptosis regulation.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.