

CHMP6 Conjugated Antibody

Catalog No: #C31838

Package Size: #C31838-AF350 100ul #C31838-AF405 100ul #C31838-AF488 100ul #C31838-AF555 100ul #C31838-AF594 100ul #C31838-AF647 100ul #C31838-AF680 100ul #C31838-AF750 100ul #C31838-Biotin 100ul #C31838-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	CHMP6 Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	WB, IF
Species Reactivity	Hu, Ms
Immunogen Description	Fusion protein of human CHMP6
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Target Name	CHMP6
Other Names	VPS20
Accession No.	Swiss-Prot#: P0DN86NCBI Protein#: BC010108
Calculated MW	23 kDa
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at -20°C/1 year

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Background

This gene encodes a member of the chromatin-modifying protein/charged multivesicular body protein family. Proteins in this family are part of the ESCRT-III (endosomal sorting complex required for transport III) which degrades surface receptors, and in biosynthesis of endosomes.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.