

KCTD17 Conjugated Antibody

Catalog No: #C31844



Package Size: #C31844-AF350 100ul #C31844-AF405 100ul #C31844-AF488 100ul #C31844-AF555 100ul #C31844-AF594 100ul #C31844-AF647 100ul #C31844-AF680 100ul #C31844-AF750 100ul #C31844-Biotin 100ul #C31844-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	KCTD17 Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IF
Species Reactivity	Hu
Immunogen Description	Fusion protein of human KCTD17
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Target Name	KCTD17
Accession No.	Swiss-Prot#: Q6PJW8NCBI Protein#: BC031038
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at -20°C/1 year

Application Details

IF: 1:50-1:200

Background

This gene encodes a protein that belongs to a conserved family of potassium channel tetramerization domain (KCTD)-containing proteins. The encoded protein functions in ciliogenesis by acting as a substrate adaptor for the cullin3-based ubiquitin-conjugating enzyme E3 ligase, and targets trichoplein, a keratin-binding protein, for degradation via polyubiquitylation. A mutation in this gene is associated with autosomal dominant myoclonic dystonia 26.?

Note: This product is for in vitro research use only and is not intended for use in humans or animals.