

IPO11 Conjugated Antibody

Catalog No: #C31857

Package Size: #C31857-AF350 100ul #C31857-AF405 100ul #C31857-AF488 100ul #C31857-AF555 100ul #C31857-AF594 100ul #C31857-AF647 100ul #C31857-AF680 100ul #C31857-AF750 100ul #C31857-Biotin 100ul #C31857-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	IPO11 Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	WB, IF
Species Reactivity	Hu, Ms
Immunogen Description	Fusion protein of human IPO11
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Target Name	IPO11
Other Names	RanBP11
Accession No.	Swiss-Prot#: Q8TEH3NCBI Protein#: BC033776
Calculated MW	113 kDa
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at -20°C/1 year

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Background

Importins, including IPO11, are a members of the karyopherin/importin-beta family of transport receptors (see KPNB1; 602738) that mediate nucleocytoplasmic transport of protein and RNA cargoes (Plafker and Macara, 2000 [PubMed 11032817]).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.